



FACILITATOR RESOURCE

RDT game

INSTRUCTIONS FOR FACILITATOR: Photocopy the page and cut into individual strips. Make enough copies to have one set of strips for every 2 trainees (e.g., if there are 10 trainees, you will need to have 5 pages/5 sets of strips). If using a blood transfer device different than capillary tube, adapt the three steps in gray on Page 1 accordingly.

Look for the expiry date at the back of the packet
Put on new gloves
Open the packet and remove the items
Write the child's name on the test tray
Open the alcohol swab and clean the finger
Open the lancet and prick child's finger to get a drop of blood
Discard the lancet in the Sharps Box
Use the capillary tube to collect the drop of blood
Use the capillary tube to put the drop of blood into the square hole A on test tray
Discard the capillary tube in the Sharps Box
Add buffer into the round hole (marked "B") on test tray.
Wait 15 minutes after adding buffer
Read/interpret test results

Dispose of all the waste in non-sharps waste containers (one for infectious waste and one for non-infectious).

Record the test results in your register